

# Potassium Bromate with Post Column Derivatization

## 1. Analytical Condition

### YL9100 HPLC System + Post Column Derivatization System

Mobile Phase : MeOH 100mL + Water 700mL + Acetic Acid 2.0g + 10% TBAH 45g

(Adjust to pH 5.0 with 10% TBAH and Water in 1L)

Flow Rate : 0.9 ml/min

Column Oven : 40°C

Detector : UV/Vis 450 nm

Injection Volume : 100uL

Column : Restek Ultra Aqueous C18 (4.6\*250mm, 5um)

### \* Post-Column Condition

System : Pinnacle PCX

Reagent : Water 700mL + HNO<sub>3</sub> 60mL + Potassium Bromide(KBr) 10g

+ ( Dissolve o-Dianisidine-2HCl 500mg in 200 ml of MeOH completely  
and produce 1L solution with water)

Reactor condition : 60 °C, 0.3mL/min

## 2. Chromatogram

