



Analysis of Residual Pesticides by HPLC

1. Analytical Condition

YL9100 HPLC

Mobile Phase : Water, ACN Gradient Program

Flow Rate : 1.5mL/min

Column Oven : 40 °C


Detector : UV/Vis 254 nm

Column : Kinetex C18 (4.6 * 100 mm, 2.6 μm)

Injection Volume : 10 μL

Gradient Program	Water	ACN
Initial	90	10
15.00	10	90
16.00	10	90
16.10	90	10
25.00	90	10

2. Chromatogram



Analysis of Residual Pesticides by HPLC

2. Chromatogram

