# ChroZen HPLC를 활용한 벤잘코늄염화물(Benzalkonium chloride) 분석

• HPLC Application

## Abstract

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### Instruments and Software

| Item               | Description  | Part No.   |
|--------------------|--|------------|
| Pump               | ChroZen HPLC Quaternary Gradient Pump with Vacuum degasser | 9421011020 |
| Autosampler        | ChroZen HPLC Autosampler                                   | 6421011020 |
| Column Compartment | ChroZen HPLC Column Compartment for Analytical scale       | 3421011020 |
| Detector           | ChroZen HPLC UV/Vis Detector with dual wavelength          | 7411011020 |
| Install. Option    | HPLC Performance Kit (Without LC C18 Column)               | 1601011890 |
| CDS                | YL-Clarity software for single instrument of YL HPLC       | 5301011000 |
|                    | Autosampler control of YL-Clarity                          | 5301011040 |
| Column             | CN (4.6 mm x 150 mm, 5 μm)                                 | -          |

### Reagents and Standards

- Acetonitrile (CH<sub>3</sub>OH), HPLC Grade
- Benzalkonium chloride, ≥ 50%
- Phosphoric acid, 85% purity
- Sodium phosphate dibasic, 99% purity
- Ultrapure water, 18.2 MΩ-cm resistivity

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


## Instrument conditions & Chromatogram

| ChroZen HPLC system |   |
|---------------------|---|
| Mobile phase        | ACN : 1% Na <sub>2</sub> HPO <sub>4</sub> = 50 : 50 (pH 5.2, phosphoric acid) |
| Column              | CN (4.6 mm x 150 mm x 5 μm)   |
| Flow rate           | 1.0 mL/min  |
| Temperature         | 40°C  |
| Injection volume    | 20 μL   |
| Detection           | UV/Vis detector 254 nm  |

| STD       | Standard Concentration (µg/mL) | Area (mV.s) |
|-----------|--------------------------------|-------------|
| 1         | 5                              | 3.5105      |
| 2         | 10                             | 7.0135      |
| 3         | 20                             | 12.8988     |
| 4         | 50                             | 36.6985     |
| 5         | 100                            | 73.6270     |
| Linearity |                                | 0.9997043   |

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### Instruments and Software

| Item               | Description  | Part No.   |
|--------------------|--|------------|
| Pump               | ChroZen HPLC Quaternary Gradient Pump with Vacuum degasser | 9421011020 |
| Autosampler        | ChroZen HPLC Autosampler                                   | 6421011020 |
| Column Compartment | ChroZen HPLC Column Compartment for Analytical scale       | 3421011020 |
| Detector           | ChroZen HPLC UV/Vis Detector with dual wavelength          | 7411011020 |
| Install. Option    | HPLC Performance Kit (Without LC C18 Column)               | 1601011890 |
| CDS                | YL-Clarity software for single instrument of YL HPLC       | 5301011000 |
|                    | Autosampler control of YL-Clarity                          | 5301011040 |
| Column             | C18 (4.6 mm x 250 mm, 5 μm)                                | -          |

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