

## 「Application Note」

# Simultaneous Analysis of DL-Camphor, L-Menthol in Pain Relief Patch by ChroZen GC

- GC Application



## Analytical Conditions

Column : DB-WAX 30 m x 0.25 mm x 0.25  $\mu$ m  
Oven : 100°C(1 min), 10°C/min, 160°C(3 min), 40°C/min,  
240°C(4 min),  
Carrier Gas : He  
Flow Rate : 2 mL/min (10:1)  
Injector : Capillary, 260°C  
Detector : FID, 260°C  
Injection Vol. : 1  $\mu$ L

## Results

